

HDF Medical Product Catalogue



Mfg. of: Maxillofacial Implants & Instruments
Orthopedic implants & Instruments

MINI SCREW SINGLE SLOT

- Product Code** : 1.5mm : TA01.04 To TA01.12
2.0mm : TA02.04 To TA02.16
2.5mm : TA03.06 To TA03.18
- Material** : Stainless Steel 316L, Titanium (Grade 5)
- Size** : 1.5mm : 4mm To 12mm
2.0mm : 4mm To 16mm
2.5mm : 6mm To 18mm



MINI SCREW CROSS SLOT

- Product Code** : 1.5mm : TA04.04 To TA04.12
2.0mm : TA05.04 To TA05.16
2.5mm : TA06.06 To TA06.18
- Material** : Stainless Steel 316L, Titanium (Grade 5)
- Size** : 1.5mm : 4mm To 12mm
2.0mm : 4mm To 16mm
2.5mm : 6mm To 18mm

INTERMAXILLARY FIXATION SCREWS (IMF) WITH GROVE & HOLE

- Product Code** : 2.0mm : TA12.04 To TA12.12
2.5mm : TA13.06 To TA13.14
- Material** : Stainless Steel 316L, Titanium (Grade 5)
- Size** : 2.0mm : 4mm To 12mm (Diff. 2mm)
2.5mm : 6mm To 14mm (Diff. 2mm)



MINI LAG SCREW (SINGLE SLOT)

- Product Code** : 2.0mm : TA07.10 to TA07.20
2.5mm : TA08.10 to TA08.20
- Material** : Stainless Steel 316L, Titanium Grade II
- Length** : 10mm to 20mm (2mm Variation)

MINI BONE PLATE WITH GAP

Product Code	: 1.5mm : TB01.02 To TB01.06 2.0mm : TC01.02 To TC01.06 2.5mm : TD01.02 To TD01.06
Screw	: 1.5mm, 2.0mm, 2.5mm
Material	: Stainless Steel 316L, Titanium (Grade II)
Thickness	: 1.5mm:-0.8mm, 2.0mm:-1.0mm, 2.5mm:-1.2mm
Holes	: 2, 3, 4, 5, 6



MINI BONE PLATE CONTINUOUS HOLE



Product Code	: 1.5mm : TB02.02 To TB02.50 2.0mm : TC02.02 To TC02.50 2.5mm : TD02.02 To TD02.50
Screw	: 1.5mm, 2.0mm, 2.5mm
Material	: Stainless Steel 316L, Titanium (Grade II)
Thickness	: 1.5mm:-0.8mm, 2.0mm:-1.0mm, 2.5mm:-1.2mm
Holes	: 2, 3, 4, 5, 6, 8, 10, 12, 16, 20, 30, 50 Holes

MINI ORBITAL PLATE WITH GAP

Product Code	: 1.5mm : TB03.04, TB03.06 2.0mm : TC03.04, TC03.06 2.5mm : TD03.04, TD03.06
Screw	: 1.5mm, 2.0mm, 2.5mm
Material	: Stainless Steel 316L, Titanium (Grade II)
Thickness	: 1.5mm:-0.8mm, 2.0mm:-1.0mm, 2.5mm:-1.2mm
Holes	: 4 Hole, 6 Hole




MINI ORBITAL PLATE CONTINUOUS HOLE



Product Code	: 1.5mm : TB04.04 To TB04.20 2.0mm : TC04.04 To TC04.20 2.5mm : TD04.04 To TD04.20
Screw	: 1.5mm, 2.0mm, 2.5mm
Material	: Stainless Steel 316L, Titanium (Grade II)
Thickness	: 1.5mm:-0.8mm, 2.0mm:-1.0mm, 2.5mm:-1.2mm
Holes	: 4, 6, 8, 10, 12, 16, 20 Holes

MINI " L " SHAPE PLATE 90' WITH GAP (LEFT & RIGHT)

- Product Code** : 1.5mm: TB05.22 To TB05.33
2.0mm: TC05.22 To TC05.33
2.5mm: TD05.22 To TD05.33
- Screw** : 1.5mm, 2.0mm, 2.5mm
- Material** : Stainless Steel 316L, Titanium (Grade II)
- Thickness** : 1.5mm:-0.8mm, 2.0mm:-1.0mm, 2.5mm:-1.2mm
- Holes** : 2x2, 2x3, 3x2, 3x3



MINI " L " SHAPE PLATE 120' WITH GAP (LEFT & RIGHT)

- Product Code** : 1.5mm: TB07.22 To TB07.33
2.0mm: TC07.22 To TC07.33
2.5mm: TD07.22 To TD07.33
- Screw** : 1.5mm, 2.0mm, 2.5mm
- Material** : Stainless Steel 316L, Titanium (Grade II)
- Thickness** : 1.5mm:-0.8mm, 2.0mm:-1.0mm, 2.5mm:-1.2mm
- Holes** : 2x2, 2x3, 3x2, 3x3



MINI " T " SHAPE PLATE WITH GAP

- Product Code** : 1.5mm: TB08.22 To TB08.33
2.0mm: TC08.22 To TC08.33
2.5mm: TD08.22 To TD08.33
- Screw** : 1.5mm, 2.0mm, 2.5mm
- Material** : Stainless Steel 316L, Titanium (Grade II)
- Thickness** : 1.5mm:-0.8mm, 2.0mm:-1.0mm, 2.5mm:-1.2mm
- Holes** : 2x2, 2x3, 3x2, 3x3



MINI BONE PLATE " I " SHAPE

- Product Code** : 1.5mm : TB10
2.0mm : TC10
2.5mm : TD10
- Screw** : 1.5mm, 2.0mm, 2.5mm
- Material** : Stainless Steel 316L, Titanium (Grade II)
- Thickness** : 1.5mm:-0.8mm, 2.0mm:-1.0mm, 2.5mm:-1.2mm
- Holes** : 6 Hole



MINI BONE PLATE " X " SHAPE

Product Code	: 1.5mm : TB11.05 2.0mm : TC11.05 2.5mm : TD11.05
Screw	: 1.5mm, 2.0mm, 2.5mm
Material	: Stainless Steel 316L, Titanium (Grade II)
Thickness	: 1.5mm:-0.8mm, 2.0mm:-1.0mm, 2.5mm:-1.2mm
Holes	: 5 Hole



MINI BONE PLATE " Y " SHAPE



Product Code	: 1.5mm : TB12.22, TB12.23 2.0mm : TC12.22, TC12.23 2.5mm : TD12.22, TD12.23
Screw	: 1.5mm, 2.0mm, 2.5mm
Material	: Stainless Steel 316L, Titanium (Grade II)
Thickness	: 1.5mm:-0.8mm, 2.0mm:-1.0mm, 2.5mm:-1.2mm
Holes	: 2 X 2, 2 X 3

MINI BONE PLATE DOUBLE " Y " SHAPE

Product Code	: 1.5mm : TB13 2.0mm : TC13 2.5mm : TD13
Screw	: 1.5mm, 2.0mm, 2.5mm
Material	: Stainless Steel 316L, Titanium (Grade II)
Thickness	: 1.5mm:-0.8mm, 2.0mm:-1.0mm, 2.5mm:-1.2mm
Holes	: 6 Hole



MINI BONE PLATE " H " SHAPE



Product Code	: 1.5mm : TB14.33, TB14.44 2.0mm : TC14.33, TC14.44 2.5mm : TD14.33, TD14.44
Screw	: 1.5mm, 2.0mm, 2.5mm
Material	: Stainless Steel 316L, Titanium (Grade II)
Thickness	: 1.5mm:-0.8mm, 2.0mm:-1.0mm, 2.5mm:-1.2mm
Holes	: 3 X 3, 4 X 4

MINI BONE PLATE " Z " SHAPE (LEFT & RIGHT)

Product Code	: 1.5mm : TB15.22L, TB15.22R 2.0mm : TC15.22L, TC15.22R 2.5mm : TD19.22L, TD19.22R
Screw	: 1.5mm, 2.0mm, 2.5mm
Material	: Stainless Steel 316L, Titanium Grade II
Thickness	: 1.5mm - 0.8mm, 2.0mm - 1mm, 2.5mm - 1.2mm
Holes	: 2X2




MINI CURVED PLATE (LEFT & RIGHT)



Product Code	: 2.0mm : TC16.04, TC16.06 2.5mm : TD23.04, TD23.06
Screw	: 2.0mm, 2.5mm
Material	: Stainless Steel 316L, Titanium (Grade II)
Thickness	: 2.0mm:-1.0mm, 2.5mm:-1.2mm
Holes	: 4 Hole, 6 Hole

MINI BONE PLATE " L " SHAPE 90° (CONTINUES HOLE) LEFT & RIGHT

Product Code	: 1.5mm : TB06.22 to TB06.33 2.0mm : TC06.22 to TC06.33 2.5mm : TD06.22 to TD06.33
Screw	: 1.5mm, 2.0mm, 2.5mm
Material	: Stainless Steel 316L, Titanium Grade II
Thickness	: 1.5mm - 0.8mm, 2.0mm - 1mm, 2.5mm - 1.2mm
Holes	: 2X2, 2X3, 3X2, 3X3



MINI ANGULAR PLATE 150° (LEFT & RIGHT)



Product Code	: 1.5mm : TB025.04, TB25.06 2.0mm : TC35.04, TC35.06 2.5mm : TD20.04, TD20.06
Screw	: 1.5mm, 2.0mm, 2.5mm
Material	: Stainless Steel 316L, Titanium (Grade II)
Thickness	: 1.5mm:-0.8mm, 2.0mm:-1.0mm, 2.5mm:-1.2mm
Holes	: 2x2, 2x3, 3x2, 3x3

MINI ORBITAL FLOOR PLATE


Product Code : TB17.S, TB17.M, TB17.L

Screw : 1.5mm

Material : Titanium Grade II

Thickness : 0.6mm

Size : Small, Medium, Large



NOSE PLATE "Y" SHAPE



Product Code : TC24

Screw : 2.0mm

Material : Titanium (Grade II)

Thickness : 1.0mm Thickness

Holes : 10 Hole

NOSE PLATE "H" SHAPE

Product Code : TC23

Screw : 2.0mm

Material : Titanium (Grade II)

Thickness : 1.0mm Thickness

Holes : 12 Hole



TRAPEZOID CONDYLE PLATE (TCP)



Product Code : TC25

Screw : 2.0mm


Material : Titanium (Grade II)

Thickness : 1.0mm Thickness

Holes : 4 Hole

3 DIMENSIONAL BONE PLATE (3D) SQUARE SHAPE

Product Code	: 1.5mm: TB26.04, TB26.06, TB26.08 2.0mm: TC36.04, TC36.06, TC36.08 2.5mm: TD21.04, TD21.06, TD21.08
Screw	: 1.5mm, 2.0mm, 2.5mm
Material	: Stainless Steel 316L, Titanium (Grade II)
Thickness	: 1.5mm:-0.8mm, 2.0mm:-1.0mm, 2.5mm:-1.2mm
Holes	: 4 Hole, 6 Hole, 8 Hole



3 DIMENSIONAL BONE PLATE (3D) RECTANGULAR SHAPE



Product Code	: 1.5mm: TB28.04, TB28.06, TB28.08 2.0mm: TC18.04, TC18.06, TC18.08 2.5mm: TD22.04, TD22.06, TD22.08
Screw	: 1.5mm, 2.0mm, 2.5mm
Material	: Stainless Steel 316L, Titanium (Grade II)
Thickness	: 1.5mm:-0.8mm, 2.0mm:-1.0mm, 2.5mm:-1.2mm
Holes	: 4 Hole, 6 Hole, 8 Hole

MICRO BURR HOLE COVER PLATE

Product Code	: 1.5mm: TB27, 2.0mm: TC26
Screw	: 1.5mm, 2.0mm
Material	: Titanium (Grade II)
Thickness	: 1.5mm: 0.8mm, 2.0mm: 1.0mm Thickness
Holes	: 5 Hole



DELTA PLATE



Product Code	: TC20
Screw	: 2.0mm
Material	: Titanium (Grade II)
Thickness	: 1.0mm Thickness
Holes	: 4 Hole

TWISTED PLATE (LEFT & RIGHT)

Product Code : TC17

Screw : 2.0mm

Material : Titanium (Grade II)

Thickness : 1.0mm Thickness

Holes : 4 Hole, 5 Hole, 6 Hole



CHIN PLATE

Product Code : TC19.22, TC19.32

Screw : 2.0mm

Material : Titanium (Grade II)

Thickness : 1.0mm Thickness

Holes : 2x2, 3x2 (10mm Bar), 2x2, 3x2 (12mm Bar)

MINI SUBCONDYLAR STRUT PLATE (TRIFIX PLATE)

Product Code : TB31.06L, TB31.06R

Screw : 2mm

Material : Titanium Grade II

Thickness : 1mm

Size : Left & Right



MINI MAXILLA PLATE



Product Code : TC33

Screw : 2.0mm


Material : Titanium Grade II

Thickness : 1.0mm

Holes : Universal

MANDIBLE WIRE MESH (STRAIGHT)

- Product Code** : 1.5mm: TB19.22 To TB19.46
2.0mm: TC27.22 To TC27.46
- Screw** : 1.5mm, 2.0mm
- Material** : Titanium (Grade II)
- Thickness** : 1.5mm:- 0.5mm, 2.0mm:- 0.6mm
- Size (In inch)** : 2x2, 2x3, 3x3, 4x4, 4x6, 5X5, 6X6, 6X8



ARCH BAR (PACKET



- Product Code** : W04
- Screw** : -----
- Material** : Stainless Steel 316L
- Thickness** : -----
- Holes** : 300 cm

WIRE REEL

- Product Code** : 265.16 to 265.30
- Dia.** : 16 To 30 SWG (Diff. 2)
- Length** : 500mm to 6000mm
(Diff. 20mm or 50mm)



MINI SELF DRILLING SCREWS (CROSS SLOT)



- Product Code** : 1.5mm : 4mm to 12mm
2.0mm : 4mm to 14mm
- Material** : Titanium Grade 5
- Size** : 4mm to 14mm

MINI PLATE WITH GAP (NEURO)

- Product Code** : TB18.02 to TB18.06
- Screw** : 1.5mm, 2mm
- Material** : Titanium Grade II
- Thickness** : 0.6mm
- Holes** : 2 Hole (10mm Bar & 12mm Bar),
3 Hole, 4 Hole, 6 Hole



MINI ORBITAL PLATE WITH GAP & CONTINUES HOLE (NEURO)



- Product Code** : With Gap : TB35.04, TB35.06
Continues Hole : TB36.04, TB36.06
- Screw** : 1.5mm, 2.0mm
- Material** : Titanium Grade II
- Thickness** : 0.6mm
- Holes** : 4 Hole, 6Hole

MINI BONE PLATE " X " SHAPE (NEURO)

- Product Code** : TB32.04
- Screw** : 1.5mm, 2mm
- Material** : Titanium Grade II
- Thickness** : 0.6mm
- Holes** : 4 Hole



MINI BONE PLATE " L " SHAPE 90° (NEURO) LEFT & RIGHT



- Product Code** : TB30.22L, TB30.22R
- Screw** : 1.5mm, 2mm
- Material** : Titanium Grade II
- Thickness** : 0.6mm
- Holes** : 2X2 Hole

MINI BONE PLATE " Y " SHAPE (NEURO)

Product Code : TB33.05

Screw : 1.5mm, 2mm

Material : Titanium Grade II

Thickness : 0.6mm

Holes : 5 Hole



MINI BONE PLATE DOUBLE " Y " SHAPE (NEURO)

Product Code : TB34.06

Screw : 1.5mm, 2.0mm

Material : Titanium Grade II

Thickness : 0.6mm

Holes : 6 Hole

MINI BURR HOLE PREBAND PLATE (NEURO)

Product Code : 1.5mm : TB24.S, TB24.L
2.0MM : TC34.S, TC34.L

Screw : 1.5mm. 2.0mm

Material : Titanium Grade II

Thickness : 0.6mm

Size : Small - 5 hole, Large - 6 Hole



MINI MASTOID PLATE (NEURO)

Product Code : 1.5mm : TB29.S, TB29.L
2.0MM : TC22.S, TC22.L

Screw : 1.5mm. 2.0mm


Material : Titanium Grade II

Thickness : 0.6mm

Size : Small & Large


EAZY MESH (Straight)

- Product Code** : 1.5mm: TB21.22 To TB21.68
2.0mm: TC29.22 To TC29.68
- Screw** : 1.5mm, 2.0mm,
- Material** : Titanium (Grade II)
- Thickness** : 0.5mm, 0.6mm, 0.8mm
- Size (In inch)** : 2x2, 2x3, 3x3, 4x4, 4x6, 5X5, 6X6, 6X8



EAZY MESH (Pre band)

- Product Code** : 1.5mm: TB22.22 To TB22.68
2.0mm: TC30.22 To TC30.68
- Screw** : 1.5mm, 2.0mm,
- Material** : Titanium (Grade II)
- Thickness** : 0.5mm, 0.6mm, 0.8mm
- Size (In inch)** : 2x2, 2x3, 3x3, 4x4, 4x6, 5X5, 6X6, 6X8



2.7MM RECONSTRUCTION PLATE (STRAIGHT)

- Product Code** : TE01.06 to TE01.50
- Screw** : 2.5mm
- Material** : Titanium (Grade II)
- Thickness** : 2.0mm
- Size (In inch)** : 6 Hole to 50 Hole



2.7MM ANGLED RECONSTRUCTION PLATE (LEFT & RIGHT)

- Product Code** : TE02.410 to TE02.424
- Screw** : 2.5mm
- Material** : Titanium Grade II
- Thickness** : 2mm
- Holes** : 4X10, 4X16, 4X20, 4X24



2.7MM DOUBLE ANGLED MANDIBLE RECONSTRUCTION PLATE

Product Code : TE03.4164 to TE03.4244

Screw : 2.5mm

Material : Titanium Grade II

Thickness : 2mm

Holes : 4X16X4, 4X20X4, 4X24X4



2.7MM TMJ CONDYLE WITH RECONSTRUCTION PLATE (LEFT & RIGHT)

Product Code : TE05.410 To TE05.424

Screw : 2.5mm

Material : Stainless Steel 316L, Titanium (Grade II)

Thickness : 2.0mm & 2.5mm

Holes : 4x10, 4x16, 4x20, 4x24



2.7MM FULL MANDIBLE RECONSTRUCTION PLATE


Product Code : TE04.12412 to TE04.14414

Screw : 2.5mm



Material : Titanium Grade II

Thickness : 2mm





Holes : 12X4X12, 14X4X14






MAXILLOFACIAL IMPORTED GRAPHIC DESIGN KIT





MAXILLOFACIAL IMPLANT ADVANCED KIT



MAXILLOFACIAL TRAUMA IMPLANT KIT




MAXILLOFACIAL INSTRUMENTS SET (CROSS SLOT)



MAXILLOFACIAL INSTRUMENTS SET (SINGLE SLOT)



IMPACTION KIT



MANDIBLE RECONSTRUCTION PLATING SET

